

BinaxNOW[®]

Influenza A & B Rapid Test

See Product Insert for complete use
instructions and performance claims

PART 1 - SAMPLE PREPARATION PROCEDURE

Nasal Wash/Aspirate:

Nasal wash/aspirate samples do not need preparation. Go to PART 2.

Nasopharyngeal and Nasal Swabs:

Elute swab in 0.5-3.0 ml of an appropriate liquid transport system by vigorously rotating the swab in the liquid. Go to PART 2. If eluting swab in the Binax elution solution, follow the Swab Elution procedure below.

Swab Elution (Control and Patient) using Binax Elution Solution:

1. The test kit contains test vials pre-filled with elution solution. Twist off the test vial cap.
2. Put swab to be tested into test vial. Rotate swab vigorously three (3) times in the liquid.
3. Press the swab against the side of the vial and turn as you remove it from the vial. This removes sample from the swab.
4. Discard the swab.
5. Test liquid sample (from the test vial) in the NOW[®] test as soon as possible. Go to PART 2.



Manufactured by:

 **inverness medical**
professional diagnostics

EC REP

EMERGO EUROPE
P.O. BOX 18510
2502 EM THE HAGUE
THE NETHERLANDS



Binax, Inc.
10 Southgate Road, Scarborough, Maine 04074
US: 1-800-323-3199 Outside US: 1-207-730-5700

BinaxNOW[®]


Influenza A & B Rapid Test



Store kits at room temperature (15°-30°C, 59°-86°F)

PART 2 - TEST PROCEDURE

1. Remove device from pouch just prior to testing and lay flat on work bench.
2. Fill pipette by firmly squeezing the top bulb and placing pipette tip into sample. Release bulb while tip is still in sample. This will pull liquid into the pipette. Make sure there are no air spaces in the lower part of the pipette.
3. See arrow on test device to find **WHITE** sample pad at the top of the test strip. **SLOWLY** (drop by drop) add entire contents of the pipette (100µl) to the **MIDDLE** of this pad such that all of the sample volume absorbs into this pad. **DO NOT** add sample to the pink/purple colored pad.
4. Immediately peel off the adhesive liner from the test device. Close and securely seal the device. Read result in window 15 minutes after closing the device.



NOTE: When reading test results, tilt the device to reduce glare on the result window if necessary.

PART 3 - RESULT INTERPRETATION

For a **NEGATIVE SAMPLE**, the BLUE Control Line in the **BOTTOM THIRD** of the window turns a pink-to-purple color. No other line appears.



For a **FLU A POSITIVE SAMPLE**, the BLUE Control Line turns a pink-to-purple color AND a second pink-to-purple Sample Line appears above it in the **MIDDLE THIRD** of the window. Any Sample Line, even when very faint, is positive.



For a **FLU B POSITIVE SAMPLE**, the BLUE Control Line turns a pink-to-purple color AND a second pink-to-purple Sample Line appears above it in the **TOP THIRD** of the window. Any Sample Line, even when very faint, is positive.



A test is **INVALID** if the Control Line remains BLUE or is not present at all, whether a Sample Line(s) is present or not. Repeat **INVALID** tests with a new test device.

